



Product Catalog

Mongoose[®] Pediatric Angiography Catheters
Super Sheath[®] Pediatric Introducers
Spring Guide Wires




FDA CLEARED FOR PEDIATRIC USE

Phone: (888) 691-6161
Fax: (888) 691-6161
E-mail: sales@PediaVascular.com

Mongoose[®]

Pediatric Angiography Catheters

FDA Cleared For Pediatric Use

Product Code	Size	Tip Shape	Length	Maximum Wire	Maximum Flow Rate*	
A-3PIG-40/0001	3.3F	Pigtail	40 cm	.030"	13 mL/sec	
A-3PIG-60/0002	3.3F	Pigtail	60 cm	.030"	11 mL/sec	
A-3PIG-80/0003	3.3F	Pigtail	80 cm	.030"	9 mL/sec	
A-4PIG-60/0004	4F	Pigtail	60 cm	.035"	22 mL/sec	
A-4PIG-80/0005	4F	Pigtail	80 cm	.035"	19 mL/sec	
A-5PIG-60/0006	5F	Pigtail	60 cm	.038"	32 mL/sec	
A-5PIG-80/0007	5F	Pigtail	80 cm	.038"	28 mL/sec	
A-3JB1-60/0008	3.3F	JB1	60 cm	.030"	11 mL/sec	
A-3JB1-80/0009	3.3F	JB1	80 cm	.030"	9 mL/sec	
A-3JB1-60/0044 ¹	3.3F	JB1	60 cm	.030"	11 mL/sec	
A-3JB1-80/0045 ¹	3.3F	JB1	80 cm	.030"	9 mL/sec	
A-4JB1-60/0010	4F	JB1	60 cm	.038"	22 mL/sec	
A-4JB1-80/0011	4F	JB1	80 cm	.038"	19 mL/sec	
A-5JB1-60/0012	5F	JB1	60 cm	.038"	32 mL/sec	
A-5JB1-80/0013	5F	JB1	80 cm	.038"	28 mL/sec	
A-3JR1-60/0014	3.3F	JR1	60 cm	.030"	11 mL/sec	
A-3JR1-80/0015	3.3F	JR1	80 cm	.030"	9 mL/sec	
A-4JR1-60/0016	4F	JR1	60 cm	.038"	22 mL/sec	
A-4JR1-80/0017	4F	JR1	80 cm	.038"	19 mL/sec	
A-5JR1-60/0018	5F	JR1	60 cm	.038"	32 mL/sec	
A-5JR1-80/0019	5F	JR1	80 cm	.038"	28 mL/sec	





¹ 0044 and 0045 JB1 catheters have no side holes, designed for placing coils

*@ 1200 PSI

Mongoose[®]

Pediatric Angiography Catheters

FDA Cleared For Pediatric Use

Product Code	Size	Tip Shape	Length	Maximum Wire	Maximum Flow Rate*	
A-3JL1-60/0020	3.3F	JL1	60 cm	.030"	11 mL/sec	
A-3JL1-80/0021	3.3F	JL1	80 cm	.030"	9 mL/sec	
A-4JL1-60/0022	4F	JL1	60 cm	.038"	22 mL/sec	
A-4JL1-80/0023	4F	JL1	80 cm	.038"	19 mL/sec	
A-5JL1-60/0024	5F	JL1	60 cm	.038"	32 mL/sec	
A-5JL1-80/0025	5F	JL1	80 cm	.038"	28 mL/sec	
A-3JLA-60/0026	3.3F	JL1.5	60 cm	.030"	11 mL/sec	
A-3JLA-80/0027	3.3F	JL1.5	80 cm	.030"	9 mL/sec	
A-4JLA-60/0028	4F	JL1.5	60 cm	.038"	22 mL/sec	
A-4JLA-80/0029	4F	JL1.5	80 cm	.038"	19 mL/sec	
A-5JLA-60/0030	5F	JL1.5	60 cm	.038"	32 mL/sec	
A-5JLA-80/0031	5F	JL1.5	80 cm	.038"	28 mL/sec	
A-3MP-60/0032	3.3F	MP	60 cm	.030"	11 mL/sec	
A-3MP-80/0033	3.3F	MP	80 cm	.030"	9 mL/sec	
A-4MP-60/0034	4F	MP	60 cm	.038"	22 mL/sec	
A-4MP-80/0035	4F	MP	80 cm	.038"	19 mL/sec	
A-5MP-60/0036	5F	MP	60 cm	.038"	32 mL/sec	
A-5MP-80/0037	5F	MP	80 cm	.038"	28 mL/sec	
A-3COB-60/0038	3.3F	Cobra	60 cm	.030"	11 mL/sec	
A-3COB-80/0039	3.3F	Cobra	80 cm	.030"	9 mL/sec	
A-4COB-60/0040	4F	Cobra	60 cm	.038"	22 mL/sec	
A-4COB-80/0041	4F	Cobra	80 cm	.038"	19 mL/sec	
A-5COB-60/0042	5F	Cobra	60 cm	.038"	32 mL/sec	
A-5COB-80/0043	5F	Cobra	80 cm	.038"	28 mL/sec	

*@ 1200 PSI

Super Sheath[®]

Pediatric Introducers

FDA Cleared For Pediatric Use



Product Code	Size	Dilator Taper	Length	Wire (kit)
S-3N5/0895	3.3F	.018"	5 cm	n/a
S-3N7/0896	3.3F	.018"	7 cm	n/a
S-3W5/0897	3.3F	.018"	5 cm	.018"
S-3W7/0898	3.3F	.018"	7 cm	.018"
S-3LONG/0899	3F	.018"	40 cm	n/a

Spring Wires

Product Code	Size	Length	Floppy Tip
SW30-140	.030"	140 cm	2 cm
SW30-140J	.030"	140 cm	5 mm J

Phone: (888) 691-6161
 Fax: (888) 691-6161
 E-mail: sales@PediaVascular.com

