


Guidewires



- **Atraumatic "J" tip design to prevent wall injury**
- **Exchange length up to 300 centimeters**
- **Softer tip and flexible shaft design allows atraumatic positioning**
- **Ultra short flexible tips**

Guidewires Offering Catheter Support and Tracking Capabilities Across a Wide Range of Vascular Anatomy

Helping Hearts Heal Themselves

AMPLATZER® Guidewires

Intended Use

AMPLATZER® Guidewires are intended for percutaneous vessel entry typically using the Seldinger technique to facilitate subsequent introduction of an intravascular device.


Description

AMPLATZER Guidewires are manufactured from stainless steel and are coated with PTFE, offering excellent glide or positioning capabilities. Guidewire length, diameter, core configuration, and tip configuration are all indicated on the product label.


Note: Refer to Instructions for Use for a complete listing of Indications and Usage, Contraindications, Warnings and Precautions. Contact AGA Medical Corporation for further information.



AMPLATZER 3 mm Modified J-Tip Guidewire in a PDA case



AMPLATZER Noodlewire in a Perimembranous VSD case




AMPLATZER 3mm Modified J-Tip Guidewire, 9-GW-001



AMPLATZER 1.5mm Modified J-Tip Guidewire, 9-GW-002



AMPLATZER 6mm Modified J-Tip Guidewire, 9-GW-003



AMPLATZER Noodlewire, 9-GW-004

AMPLATZER Guidewire Ordering Information

Order No.	Wire Size	Type	Tip Description	Usable Length
9-GW-001	.035"	Super Stiff	3 mm Modified J-Tip	260 cm
9-GW-002	.035"	Super Stiff	1.5 mm Modified J-Tip	260 cm
9-GW-003	.035"	Super Stiff	6 mm Modified J-Tip	300 cm
9-GW-004	.035"	Soft Fixed Core	6 mm J-Tip	300 cm

The AMPLATZER Guidewires are Designed to Work With the Following AMPLATZER Occluders

9-GW-001 Guidewire	AMPLATZER PDA Occluder
9-GW-002 Guidewire	AMPLATZER ASD Occluder
9-GW-003 Guidewire	AMPLATZER Muscular VSD Occluder
9-GW-004 Noodlewire	AMPLATZER Membranous VSD Occluder



AGA Medical Corporation
682 Mendelssohn Avenue
Golden Valley, MN 55427 USA

(888) 546-4407 Toll Free
(763) 513-9227 Phone
(763) 513-9226 Fax
www.amplatzer.com Internet

AMPLATZER® is a registered trademark of AGA Medical Corporation. TorqVue™ is a trademark of AGA Medical Corporation. ©AGA MEDICAL CORPORATION, 2003. Manufactured under one or more of the following patents: U.S. Patent 5,725,552; 5,846,261; 5,944,738; 6,123,715 & Other U.S. and Foreign patents pending. Not in any way connected with medical gas or equipment sold under the "AGA" brand by AGA AB or its successors.